



## Laser Diode Controller

# LDC200C Series Operation Manual



2018

---

Version: 6.5  
Date: 10-Jul-2018

# Contents

<b>Foreword</b>	<b>3</b>
<b>1 General Information</b>	<b>4</b>
1.1 Laser Diode Protection Features	5
1.2 Safety	5
1.3 Ordering Codes and Accessories	7
<b>2 Getting Started</b>	<b>8</b>
2.1 Parts List	8
2.2 Preparation	8
2.3 Operating elements	9
2.4 First Operation	10
<b>3 Operating Instruction</b>	<b>11</b>
3.1 External Connections	11
3.1.1 Laser Diode Output	11
3.1.2 LD and PD polarities	12
3.1.3 Interlock	13
3.1.4 Analog Out	14
3.2 Setting the Current Limit	14
3.3 Constant Current Mode (CC)	14
3.4 Constant Power Mode (CP)	15
3.5 Calibration of Optical Power Display	16
3.6 External Analog Modulation	17
3.7 Overtemperature Protection	17
3.8 Disabling the Beeper	18
3.9 Remote Operation	18
<b>4 Maintenance and Service</b>	<b>20</b>
4.1 Line Voltage Setting	20
4.2 Replacing Mains Fuses	21
4.3 Troubleshooting	22
<b>5 Appendix</b>	<b>24</b>
5.1 Technical Data	24
5.1.1 LDC200CV	25
5.1.2 LDC201CU	26
5.1.3 LDC202C	27
5.1.4 LDC205C	28
5.1.5 LDC210C	29
5.1.6 LDC220C	30
5.1.7 LDC240C	31
5.2 Declaration of Conformity	32

<b>5.3 Warranty</b>	<b>33</b>
<b>5.4 Copyright and Exclusion of Reliability</b>	<b>34</b>
<b>5.5 Thorlabs 'End of Life' Policy (WEEE)</b>	<b>35</b>
<b>5.6 Thorlabs Worldwide Contacts</b>	<b>36</b>

We aim to develop and produce the best solution for your application in the field of optical measurement technique. To help us to live up to your expectations and improve our products permanently we need your ideas and suggestions. Therefore, please let us know about possible criticism or ideas. We and our international partners are looking forward to hearing from you.

*Thorlabs GmbH*

### **Warning**

Sections marked by this symbol explain dangers that might result in personal injury or death. Always read the associated information carefully, before performing the indicated procedure.

### **Attention**

Paragraphs preceded by this symbol explain hazards that could damage the instrument and the connected equipment or may cause loss of data.

### **Note**

This manual also contains "NOTES" and "HINTS" written in this form.

Please read these advices carefully!

# 1 General Information

The Thorlabs LDC200C Series Laser Diode Controllers are high accuracy precise injection current controllers for laser diodes and LEDs. Together with a Thorlabs Temperature Controller a stable operation of the connected laser diode can be achieved. The LDC200C Series includes the following types:

- LDC200CV - designed for safe operation of VCSEL laser diodes.
- LDC201CU - ultra low noise current ( $<0.2\mu\text{A}$  RMS).
- LDC202C, LDC205C and LDC210C - enhanced compliance voltage ( $>10\text{V}$ ) for use with blue laser diodes.
- LDC240C - higher current (4A).

The LDC200C Series controllers are easy to operate via the operating elements on the front panel. Operating parameters are shown on a 5-digit LED display. UP-DOWN keys allow to select the parameter to be displayed.

After switching on a LDC200C Series laser diode controller, it remains in LASER OFF mode. The laser current can be switched on/off using the appropriate key at the front panel.

Additionally the laser current can be switched by applying a TTL signal to the LD remote input at the rear of the unit.

The laser and the photodiode are connected via a 9-pin D-SUB jack at the rear of the unit. The output for the laser diode and the input for the photodiode are bipolar, thus all polarities of commercial available laser diodes can be connected.

The injection current or the optical output power of the laser diode can be modulated applying a modulation signal to the input at the rear of the unit.

A voltage proportional to the laser diode current is provided for monitoring purposes at an analog control output at the rear.

If an error occurs or the limit for the laser current is reached, the corresponding LED lights up and a short beep gives a warning.

For a low ripple and noise of the output current a mains filter is installed and the transformer is shielded carefully.

The LDC200C Series controller are cooled by an internal fan, which protects the unit against overheating in case of high environmental temperatures. With free air circulation a safe operation of the unit is guaranteed up to 40 °C ambient temperature.

## **Warning**

**Do not obstruct the air ventilation slots in the housing!**

## **Note**

**In order to prevent damages to the laser diode, it is recommended to mount the laser into a suitable Thorlabs laser diode mount and connect it to the LDC200C Series using the supplied Thorlabs CAB400 cable. This ensures the utmost protection of the laser diode from damage by wrong connection.**

## 1.1 Laser Diode Protection Features

- **Soft start function**  
The soft start function protects the laser diode against undesired transients by smoothly increasing the injection current after switching the laser on. This current ramp last for abt. 200 msec.
- **Adjustable limitation of laser current**  
An adjustable hardware limit for the laser current can be set using the 12 turn potentiometer at the front panel. This protects the laser diode from unintentional laser current values.
- **Contact protection of the laser diode (open circuit)**  
If the connection to the laser is interrupted during operation even for a very short time, the output is switched off immediately and can be switched on manually only.
- **Electronic short-circuit for the laser diode in LASER OFF mode**  
With the output switched off the laser diode is short-circuited electronically, thus the laser diode is protected against accidental ESD and induced current transients.
- **Line failure protection**  
The laser current is switched off immediately if a power failure or line interruption occurs. In this case, same as after turning on the controller, the laser current remains switched off and can be switched on manually only.
- **Interlock**  
The interlock interface provides a number of several protection functions simultaneously.
  - Safety lock to prevent unintentional use
  - Cable disconnect or damage monitoring
  - An external emergency switch can be connected
  - Application of external automatic protection and alerting equipment
  - External LED to indicate "Laser ON"

The laser can only be switched on with the interlock input closed. If the interlock contact is opened during operation (even for a very short time) the output is switched off immediately and can be switched on manually only.

## 1.2 Safety

### **Attention**

**The safety of any system incorporating the equipment is the responsibility of the assembler of the system.**

**All statements regarding safety of operation and technical data in this instruction manual will only apply when the unit is operated correctly as it was designed for.**

**The LDC200C Series must not be operated in explosion endangered environments!**

**Do not obstruct the air ventilation slots in the housing!**

**Do not remove covers!**

**Do not open the cabinet. There are no parts serviceable by the operator inside!**

**This precision device is only serviceable if properly packed into the complete original packaging including the plastic foam sleeves. If necessary, ask for replacement packaging.**

**Refer servicing to qualified personnel!**

**Only with written consent from Thorlabs may changes to single components be made or components not supplied by Thorlabs be used.**

**Attention**

Prior to applying power to the LDC200C Series, make sure that the protective conductor of the 3 conductor mains power cord is correctly connected to the protective earth ground contact of the socket outlet! Improper grounding can cause electric shock resulting in damage to your health or even death!

Ensure that the line voltage setting of the fuse holder at the rear panel agrees with your local supply and that the corresponding fuses are inserted. If not, please change the line voltage setting (see section [Line voltage setting](#)) and the mains fuses (see section [Replacing the main fuses](#)).

To avoid risk of fire, only the appropriate fuses for the corresponding line voltage must be used.

All modules must only be operated with proper shielded connection cables.

**Attention**

The following statement applies to the products covered in this manual, unless otherwise specified herein. The statement for other products will appear in the respective accompanying documentation.

**Note:**

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and meets all requirements of the Canadian Interference-Causing Equipment Standard ICES-003 for digital apparatus. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/T.V. technician for help.

Users that change or modify the product described in this manual in a way not expressly approved by Thorlabs (party responsible for compliance) could void the user's authority to operate the equipment.

Thorlabs GmbH is not responsible for any radio television interference caused by modifications of this equipment or the substitution or attachment of connecting cables and equipment other than those specified by Thorlabs GmbH. The correction of interference caused by such unauthorized modification, substitution or attachment will be the responsibility of the user.

The use of shielded I/O cables is required when connecting this equipment to any and all optional peripheral or host devices. Failure to do so may violate FCC and ICES rules.

**Attention**

Mobile telephones, cellular phones or other radio transmitters must not be used within the range of three meters of this unit since the electromagnetic field intensity may then exceed the maximum allowed disturbance values according to IEC61326-1.



This product has been tested and found complying with the limits according to IEC 61326-1 for using connection cables shorter than or equal to 3 meters (9.8 feet).

### Warning

Laser modules can deliver up to several 100mW of even invisible laser radiation! When operated incorrectly, this can cause severe damage to your eyes and health! Be ensure to pay strict attention to the safety recommendations of the appropriate laser safety class, as stated for the used light source.

Vous pouvez trouver les traductions françaises des paragraphes ayant trait à la sécurité d'utilisation de ce produit sur le lien suivant:

<https://www.thorlabs.com/sd.cfm?fileName=15988-D03.pdf&partNumber=LDC202C>

En outre, vous pouvez soit scanner le QR code, soit vous référer à la section "Documents" sur la page web du produit.



## 1.3 Ordering Codes and Accessories

### Ordering code Short description

<b>LDC200CV</b>	Laser diode controller, current range 0 ... 20 mA / 6 V
<b>LDC201CU</b>	Laser diode controller, current range 0 ... 100 mA / 5 V, Ultra Low Noise
<b>LDC202C</b>	Laser diode controller, current range 0 ... 200 mA / 10 V
<b>LDC205C</b>	Laser diode controller, current range 0 ... 500 mA / 10 V
<b>LDC210C</b>	Laser diode controller, current range 0 ... 1 A / 10 V
<b>LDC220C</b>	Laser diode controller, current range 0 ... 2 A / 4 V
<b>LDC240C</b>	Laser diode controller, current range 0 ... 4 A / 5 V

### Shielded cable:

<b>CAB400</b>	Cable to connect the LDC200C Series controller to a Thorlabs Laser Diode Mount (included)
---------------	---

### Laser diode mounts for different laser diode packages:

<b>TCLDM9</b>	Temperature controlled laser diode mount for 3- and 4-pin TO18-packages (9 mm CD, 5.6 mm CD)
<b>LDM21</b>	Miniature sized temperature controlled laser diode mount for 3- and 4-pin TO18-packages (9 mm CD, 5.6 mm CD)
<b>LM14S2</b>	laser diode mount for laser modules in a 14-pin butterfly-package (programmable pinning)

Please visit our homepage <http://www.thorlabs.com> for further information.

## 2 Getting Started

### 2.1 Parts List

Inspect the shipping container for damage.

If the shipping container seems to be damaged, keep it until you have inspected the contents and you have inspected the LDC200C Series controller mechanically and electrically.

Verify that you have received the following items within the package:

1. 1 LDC200C Series controller
2. 1 power cord, connector according to ordering country
3. 1 operation manual
4. 1 connection cable CAB400

### 2.2 Preparation


Prior to operate a LDC200C Series controller, check if the set line voltage matches with your local power supply and if the appropriate fuses are inserted. (See sections [Line Voltage Setting](#) and [Replacing the Mains Fuses](#) )

Connect the unit to the power line using the supplied cable. Turn the unit on by pressing the [line switch \(F10\)](#).

If required, the chassis ground can be connected to ground potential via the connector jack (R5). The ground pin of the laser diode is internally connected to chassis ground.


## 2.3 Operating elements

### Front Panel



<b>F1</b>	-	5-digit LED display
<b>F2</b>	LED "mA"	Current display in mA
<b>F3</b>	LED "mW"	Power display in mW
<b>F4</b>	LED "OPEN"	No laser diode connected, or Interlock open
<b>F5</b>	LED "LIMIT"	Adjusted current limit reached
<b>F6</b>	LED "OTP"	Overtemperature protection is active
<b>F7</b>	LED "LASER ON"	Laser current is switched on
<b>F8</b>	Key "LASER ON"	On / Off switch for the laser current
<b>F9</b>	-	Knob for adjusting the current or power set value
<b>F10</b>	-	Line switch (ON / OFF)
<b>F11</b>	LIM I	Potentiometer for setting the current limit
<b>F12</b>	CAL	Potentiometer for calibrating the power display
<b>F13</b>	LED "ILD"	Display shows the laser current
<b>F14</b>	LED "PLD"	Display shows the optical power
<b>F15</b>	LED "IPD"	Display shows the photodiode current
<b>F16</b>	LED "ILIM"	Display shows the current limit
<b>F17</b>	Key "DOWN"	Select the parameter to be displayed
<b>F18</b>	Key "UP"	Select the parameter to be displayed
<b>F19</b>	LED "AG"	Selected laser polarity: anode grounded
<b>F20</b>	LED "CG"	Selected laser polarity: cathode grounded
<b>F21</b>	LED "P"	Constant power mode
<b>F22</b>	LED "I"	Constant current mode
<b>F23</b>	Key "LD POL"	Select laser polarity: anode grounded or cathode grounded
<b>F24</b>	Key "CONST"	Select constant current mode or constant power mode
<b>F25</b>	PD RANGE	Potentiometer for setting the photodiode current range

## Rear Panel



- R1 TTL input "LD REM", 0 ... +5 V
- R2 Modulation input / analog control input "MOD IN", -10V ... +10 V
- R3 Analog monitoring output "CTL OUT", 0 ... ±10V
- R4 Fan
- R5 4 mm banana jack for chassis ground
- R6 Connector "LD OUT" for laser diode, photodiode, interlock, status LED
- R7 Serial number of the unit
- R8 Indicator / switch for line voltage (included in fuse holder)
- R9 Mains connector and fuse holder

## 2.4 First Operation

### Attention

Prior to switching on your LDC200C Series controller, please make sure that the [line voltage setting](#) corresponds to your mains voltage! Mismatching may lead to damage of the controller!

Turn the unit on by pressing the line switch (F10).

After switching on the unit, the LED display (F1) and a LED, indicating the selected measurement value (F13 ... F16), light up, otherwise please check the [line voltage](#) and the [mains fuses](#).

By using the keys (F17) and (F18) the displayed measurement value can be selected.

A LDC200C Series controller is immediately ready to use after turning on. The rated accuracy is reached, however, after a warming-up time of approx. 10 minutes.

## 3 Operating Instruction

### 3.1 External Connections

Prior to switch on the LDC200C Series controller, all required external connections must be made properly. Please read the following sections carefully.


#### 3.1.1 Laser Diode Output

The Thorlabs LDC200C Series controller can drive all laser diodes depending on the model up to a maximum current of 4 A.

If a Thorlabs laser mount is used, just connect the 9-pin D-Sub jack "LD OUT" (R6) of LDC200C Series controller to the 9-pin plug "LD DRIVER" of the Laser Diode Mount using the supplied shielded cable CAB400.

If a Thorlabs laser mount TCLDM9 or LDM21 is used, in addition the polarity of laser diode and photodiode must be set using the two slide switches at the laser mount. Please refer to the individual operation manual of the laser diode mount.

If other laser diode fixtures are used, connect the laser diode and - if provided - the photodiode using shielded cables to "LD OUT" (R6) according to the pin assignment as shown below:



*Pin assignment of the "LD OUT" jack (female, rear panel view)*

#### Pin # Connection

##### Interlock (open circuit monitoring, status display)

- |   |                                     |
|---|-------------------------------------|
| 1 | interlock and status LASER ON / OFF |
| 5 | ground for pin 1                    |

##### Laser Diode

- |   |   |
|---|---|
| 7 | laser diode cathode (in case of polarity "anode grounded" - AG) |
| 8 | laser diode anode (in case of polarity "cathode grounded" - CG) |
| 3 | laser diode ground  |

##### Photodiode

- |   |                    |
|---|--------------------|
| 2 | photodiode cathode |
| 4 | photodiode anode   |

#### Note

Do not connect the laser diode between pin 7 and 8! The correct connection must be:


- for laser diode "AG": cathode to pin 7, anode to pin 3, pin 8 - open
- for laser diode "CG": anode to pin 8, cathode to pin 3, pin 7 - open.

### 3.1.2 LD and PD polarities

Connect laser and photodiode to the LDC200C Series controller as shown below.

The ground connection of the laser diode (pin 3) can be connected to the anode of the photodiode (pin 4), or to the cathode of the photodiode (pin 2). This connection should be as close as possible to the laser diode to avoid measurement errors.

If the selected polarity (key "LD POL" [F23](#)) and the pin connection of the laser diode do not match, no current is flowing through the laser diode.



The laser diode is always sourced with respect to ground. Compared to a floating ground scheme, the advantages of this operation mode are higher safety for the laser diode, better stability of the laser current and less introduced interferences.

Prior to turning on the laser current, the correct polarity of the laser diode must be selected:

- Make sure the laser is switched off (key "LASER ON" F8).
- Select the required polarity by pressing the key "LD POL" (F23).
- For laser diodes with a built-in monitor diode the common pin of laser diode and photodiode is the ground pin.
- The selected polarity of the output is indicated by the LED "AG" (F19) or "CG" (F20).

If the laser diode was connected correctly (see above) but the polarity of the laser diode was selected wrong with the key "LD POL" (F23), an "OPEN" error will be indicated (LED "OPEN" F4).


After pressing "LASER ON" (F8) the LED "OPEN" (F4) extinguishes and the output can be switched on by pressing "LASER ON" (F8) again.

#### Biasing the photodiode

If the photodiode shall be operated with bias voltage, a battery can be connected in series to the photodiode.

If the photodiode is of AG polarity, connect the positive terminal of the battery to the cathode of the photodiode, and the negative terminal of the battery to pin 2.

If the photodiode is of CG polarity, connect the negative terminal of the battery to the anode of the photodiode, and the positive terminal of the battery to pin 4.

**Attention**

A wrong polarity of the battery may destroy the photodiode due to a forward current flowing through it!

**3.1.3 Interlock**

Pin 1 and pin 5 of the ["LD OUT" jack \(R6\)](#) are the interface of the interlock circuit.



The interlock output represents a current source, while the voltage across the external circuit, connected to pin 1 and 5, is being observed. As soon as this voltage rises above a certain threshold ( $\sim 2.5\text{V}$ ), the external circuit is considered as "open", the laser current output is disabled and the [LED "OPEN" \(F4\)](#) lights up .

Pin 1 and pin 5 must be connected externally by a short cut wire, by a circuit (total resistance  $< 430\ \Omega$ ) or by a LED (anode to pin 1, cathode to pin 5, connect in parallel resistor  $R=1\ \text{k}\Omega$ ). A connected LED lights up when the laser current is switched on (LASER ON).

**Note**

Do not use blue LED due to their high forward voltage

Thorlabs Laser Diode Mounts use the interlock current for indication "LD ON". When using a Thorlabs TCLDM9 or LM14S2 laser diode mount, interlock functionality can be easily used for emergency switch-off of the laser by connecting an opening contact to the jack located on the laser mount. For details, please refer to the appropriate laser diode mount manual.

### 3.1.4 Analog Out

The output "CTL OUT" ([R3](#)) delivers a DC voltage, proportional to the actual laser current. The output voltage ranges from:

- 0 ...+10 V with cathode grounded
- 0 ... -10 V with anode grounded.

When connecting a load avoid ground loops. The load resistance should be >10 kΩ. The conversion coefficient for this output can be calculated to:

$$k_{CTL\ OUT} = 10V/I_{max} [V/mA]$$

where  $I_{max}$  is the max. output current for the given controller.

#### Example:

For an LDC202C with a max. output current of 200mA the conversion coefficient  $k_{CTL\ OUT}$  is equal to 0.05 V/mA.

## 3.2 Setting the Current Limit

Prior to switch on the laser diode, always set an appropriate injection current limit to protect the laser diode from destruction:

- Switch on the unit - button "LINE" (F10).
- Make sure the laser current is switched off
- Select parameter "LIM" (key F17 or F18).
- Use a screwdriver to set the desired current with the potentiometer "LIM I" (F11).

#### Note

The current limit can be displayed at any time by selecting the parameter "LIM"

If the laser current reaches it's limit during operation, the LED "LIMIT" (F5) lights up and a short beep alerts about that. The current limit can also be adjusted when the output is switched on. Avoid any adjustments when the LED "LIMIT" (F5) lights up.

## 3.3 Constant Current Mode (CC)

- Switch on the LDC2xxC.
- Make sure, that the current limit "LIM" is [adjusted properly](#).
- Make sure that the correct [laser diode polarity](#) is set.
- [Connect the laser diode](#).
- Select the display "LD" using the keys (F17) or (F18).
- Turn the adjust knob (F9) completely counter clockwise.
- Select constant current mode ("I") by pressing the key "CONST" (F24) until LED "I" (F22) lights up.
- Switch on the output by pressing the key "LASER ON" (F8). The LED "LASER ON" (F7) lights up, the output is activated and the current slowly rises (about 1 s) to the value set with knob (F9).

#### Note

The current output can be switched on only if the jack "LD OUT" (R6) is connected correctly  
The display shows the injection current "LD".



With the adjust knob (F8) the laser current can be set continuously between 0 mA and the selected current limit "I<sub>LIM</sub>".

- ⇒ If the injection current "I<sub>LD</sub>" reaches the set current limit "I<sub>LIM</sub>" during operation, the LED "LIMIT" (F5) lights up, a short alert beep is heard, and the laser current is limited to the value of the current limit "I<sub>LIM</sub>". In this case, noise and ripple do not longer correspond to the specifications for normal operation. However exceeding the set current limit "I<sub>LIM</sub>" is impossible.
- ⇒ If the connection to the laser diode is interrupted during operation, the laser diode current is switched off automatically. LED "LASER ON" (F7) extinguishes, LED "OPEN" (F4) lights up and a short beep is heard.
- ⇒ If the output is switched on while the interlock is closed and there is no laser diode connected, at first the LED "LASER ON" (F7) lights up and the output is switched on. Then the LDC2xxC recognizes the missing laser diode and switches the output off. The LED "LASER ON" (F7) extinguishes and the LED "OPEN" (F4) lights up. Connect the laser and press "LASER ON" (F8) - the LED "OPEN" (F4) extinguishes and the output switches on.
- ⇒ If a photodiode is connected, the display can be set to show the photodiode current "I<sub>PD</sub>" or the optical power "PLD" by pressing the key (F17) or (F18).

The laser current can be [modulated](#) via the connector "MOD IN" (R2).

### 3.4 Constant Power Mode (CP)

If the laser includes a photodiode (also known as monitor diode), the laser can be operated also in constant power mode. In contrast to the constant current mode, here the photodiode current is kept constant by controlling the laser injection current.

- Start operating the LDC2xxC in constant current mode (refer to chapter 3.4)
- Select the display "I<sub>PD</sub>" with the key (F17) or (F18).
- Check whether an appropriate photodiode current "I<sub>PD</sub>" is available. A photodiode current of at least the minimum control current must be available for a stable operation of the power control:
  - 5 µA for LDC200CV,
  - 25 µA for LDC201CU, LDC202C, LDC205C,
  - 50 µA for LDC210C, LDC220C, LDC240C.

#### Note

If the sign of the photodiode current "I<sub>PD</sub>" is negative, the photodiode must be reversed.

- Switch off the laser current by pressing the key "LASER ON" (F8) .
  - Select constant power mode ("P") by pressing the key "CONST" (F24) until LED "P" (F21) lights up.
  - Set the adjust knob (F9) completely counter clockwise.
  - Switch on the output by pressing "LASER ON" (F8). The LED "LASER ON" (F7) lights up, the output is activated and the current slowly increases (about 1 s) to the set value.
  - The photodiode current "I<sub>PD</sub>" can be set using (F9) and with that the optical power of the laser diode increases until the laser current "I<sub>LD</sub>" reaches the selected current limit "I<sub>LIM</sub>".
- ⇒ If the desired photodiode current "I<sub>PD</sub>" cannot be set with a sufficient accuracy (F9), the range of the knob (F9) can be fitted to the current range of the connected photodiode using the potentiometer "PD RANGE" (F25) - see below.
  - ⇒ If no photodiode is connected or the polarity of the photodiode is set wrong, the laser current "I<sub>LD</sub>" increases to the limit "I<sub>LIM</sub>" after the output is switched on.

- ⇒ If the connection to the photodiode is interrupted in constant power mode, the laser current "I<sub>LD</sub>" increases to the set current limit "I<sub>LIM</sub>".
- ⇒ If the injection current "I<sub>LD</sub>" is limited by the selected current limit "I<sub>LIM</sub>", ripple and noise do not longer correspond to the specifications for normal operation. However, the set maximum current "I<sub>LIM</sub>" cannot be exceeded.
- ⇒ If the connection to the laser diode is interrupted during operation, the output is switched off automatically. LED "LASER ON" (F7) extinguishes, LED "OPEN" (F4) lights up and a short beep is heard.
- ⇒ If the output is switched on while the interlock is closed and there is no laser diode connected, at first the LED "LASER ON" (F7) lights up and the output switches on. Then the LDC2xxC recognizes the missing laser diode, switches the output off, the LED "LASER ON" (F7) extinguishes and the LED "OPEN" (F4) lights up. Connect the laser, then press "LASER ON" (F8) - the LED "OPEN" (F4) extinguishes and the laser diode current output switches.
- ⇒ The laser current can be modulated via the connector "MOD IN" (R3) - see section [External Analog Modulation](#).

### Changing the "I<sub>PD</sub>" setting range

With potentiometer "PD RANGE" (F25) the full scale of the adjustment knob (F9) can be set in order to ease the setting of the operating point in constant power mode even with low photodiode current. The photodiode current range depends on the type of the LDC2xxC controller:

- LDC200CV: 0.6 mA to 2 mA
  - LDC201CU, LDC202C, LDC205C: 3 mA to 10 mA
  - LDC210C, LDC220C, LDC240C: 6 mA to 20 mA
- Turn the potentiometer "PD RANGE" (F25) completely counter clockwise.
  - Turn the adjust knob (F9) completely counter clockwise.
  - Select constant power mode "P" (key "CONST" F24), connect the laser diode and the photodiode and switch on the laser current by pressing the key "LASER ON" (F8).
  - Set the desired photodiode current with the knob (F9).

If the desired operating point cannot be reached even if the knob (F9) is turned completely clockwise (i.e.  $I_{PD} > 0.6$  mA), turn the potentiometer "PD RANGE" (F25) clockwise until the desired setting range has been reached.

#### Note

A change of the "PD RANGE" (F25) leads to a change of the power set value in CP mode. In constant current mode the potentiometer "PD RANGE" is without function.

## 3.5 Calibration of Optical Power Display

Additionally to the photodiode current "I<sub>PD</sub>", the optical power of the laser diode can be displayed. To achieve a correct optical power display, a calibration (potentiometer "CAL" - F12) must be carried out:

- ⇒ Select constant current mode "I" - key "CONST" (F24).
- ⇒ With the knob (F8), set the laser current "I<sub>LD</sub>" to an operating point for which the optical output power of the laser diode is known, e.g. from the data sheet of the laser diode. You may also find this operating point by displaying the monitor current "I<sub>PD</sub>" (F15). A more accurate way is to measure the laser power using an optical power meter.

- ⇒ Select the display "PLD" (F14) - keys (F17) or (F18).
- ⇒ Turn the potentiometer "CAL" (F12) completely counter clockwise.
- ⇒ Select the range at which the displayed power is least higher than the actual optical power:
- ⇒ Press the key "CONST" (F24). While holding it, press the key "Down" (F17). With every key-stroke, the decimal point of the display moves to the left. Pressing the key "Up" (F18), the decimal point moves to the right. After reaching the highest range, the decimal point jumps back to the left-most position.
- ⇒ Calibrate the display "PLD" to the value of the actual optical power by adjusting the potentiometer "CAL" (F12) clockwise.

### 3.6 External Analog Modulation

To generate a time dependent injection current "ILD" or photodiode current "IPD", these settings can be modulated via an independent ground-symmetric modulation input "MOD IN" (R2). Maximum allowed input voltage is  $-10\text{ V} \dots +10\text{ V}$ , input resistance is  $>10\text{ k}\Omega$ .

"ILD" and "IPD" are calculated as:

$$\text{ILD} = \text{ILD SET} + \text{ILD MAX} * \text{UMOD} / 10\text{ V} \quad (\text{in constant current mode})$$

$$\text{or } \text{IPD} = \text{IPD SET} + \text{IPD MAX} * \text{UMOD} / 10\text{ V} \quad (\text{in constant power mode})$$

with:

ILD MAX:	Maximum laser current (e.g. 500 mA for LDC205C)
IPD MAX:	Max. photodiode current (e.g. 10 mA for LDC205C; see Control Ranges in " <a href="#">Technical Data</a> ")
ILD SET or IPD SET:	value set with knob (F9)
UMOD:	voltage at input "MOD IN" (R2)

Start operation in constant current or constant power mode (refer to chapter 3.4 or 3.5) and adjust the desired set value with the knob (F9).

Connect the modulation source to the jack "MOD IN" (R2). Avoid ground loops when connecting the function generator.

The laser diode current "ILD" can be monitored at the analog output "[CTL OUT](#)" (R3) .

If the injection current "ILD" reaches the current limit "ILIM" in operation, the LED "LIMIT" (F5) lights up, a short beep is heard, and the laser current is limited to the value of the current limit "ILIM". In this case, ripple and noise do no longer correspond to the specifications for normal operation. However the set maximum current "ILIM" cannot be exceeded.

### 3.7 Overtemperature Protection

The LDC200C Series controllers come with an internal over-temperature protection. If the internal heat sink is overheated, the output of the controller is disabled automatically. The LED "OTP" (F6) lights up and a short beep is heard. The laser diode current is switched off immediately. Pressing the key "LASER ON" (F8) has no effect in this case.

After the internal heat sink's temperature decreased for about  $10^{\circ}\text{C}$ , the LED "OTP" (F6) extinguishes and the laser current output can be switched on again.

## 3.8 Disabling the Beeper

If audible signals are unwanted, the beeper can be disabled in this way:

- Press and hold the key “UP” (F18).
- Press the key “Down” (F17). Now the beeper state is displayed:
  - “Sd.On” - Sound ON
  - “Sd.OFF” - Sound OFF

To change the beeper state, hold the key “UP” pressed and toggle the beeper state by pressing “DOWN” key.

## 3.9 Remote Operation

Although the LDC200C Series Controllers do not provide a standard computer interface, a limited remote operation is possible using the TTL input or the analog inputs and outputs.

The following control or read out functions are available:

### LD REM

This digital input allows to switch on/off the laser independent of the key “LASER ON”. A rising edge at this input switches the laser on (if being in off state), a falling edge switches the laser off (if being in on state). The threshold voltages are TTL compatible and feature hysteresis.

When using the “LD REM” input, the key “LASER ON” remains fully functional.

### Interlock

If a relay contact or an open collector transistor is inserted into the interlock line (refer to section "[Interlock](#)"), the output can be switched off remotely at any time.

### Note

The interlock is usable only for switching off. After the interlock path was interrupted, the LDC200C Series Controller can be switched on again only manually, provided the Interlock is closed.

### Setting the laser diode current

The analog modulation input “MOD IN” input (R2) is DC coupled; the input voltage ranges from -10V to +10V. The input resistance of the input “MOD IN” is >10 kΩ.

A voltage applied to this input is converted to a current value, which is then superposed to the LD set value, adjusted with the rotary knob F9. The conversion coefficient depends on the max. nominal output current of the controller and can be calculated to

$$k_{MOD} = I_{max} / 10V [mA/V]$$

where  $I_{max}$  is the max. output current for the given controller.

Example:

For a LDC202C with a max. output current of 200mA the conversion coefficient  $k_{MOD}$  is equal to 20mA/V.

If the adjusting knob (F9) is set to it's left stop position, corresponding to a zero injection current, the laser diode set current value can be controlled externally by an analog DC voltage of (0...+10) V applied to the “MOD IN” input (R2).

If knob F9 is at its right stop position (corresponding to the max. injection current), a DC voltage of (-10...0) V can be used to decrease the injection current.

**Note**

The laser current output must be switched on, otherwise a current control via MOD IN is impossible.

**Monitoring the laser current**

The output "CTL OUT" (R3) delivers a DC voltage, proportional to the actual laser current. The output voltage ranges from:

- 0 ... +10 V with cathode grounded
- 0 ... -10 V with anode grounded.

Load resistance should be >10 kΩ. The conversion coefficient for this output is the multiplicative inverse of  $k_{MOD}$ :

$$k_{CTL\ OUT} = 10V/I_{max} \text{ [V/mA]}$$

**Note**

All operating elements of the LDC2xxC are active at any time. For remote operation, make sure that the manual settings are not changed during operation.

## 4 Maintenance and Service

Protect the LDC200C Series from adverse weather conditions. The LDC200C Series is not water resistant.

### Attention

**To avoid damage to the instrument, do not expose it to spray, liquids or solvents!**

The unit does not need a regular maintenance by the user. If necessary the unit and the display can be cleaned with a cloth dampened with water. A mild 75% Isopropyl Alcohol solution can be used for more efficient cleaning.

LDC200C Series Controllers do not contain any modules and/or components that could be repaired by the user himself. If a malfunction occurs, please contact Thorlabs for return instructions.

Do not remove covers!

To guarantee the specifications given in section [Technical Data](#) over a long period it is recommended to have the unit factory calibrated every two years.

### 4.1 Line Voltage Setting

The laser diode controller LDC2xxC operates at fixed line voltages of


100 V +15% / -10% ( 90 V ... 115 V)

115 V +15% / -10% (104 V ... 132 V)

230 V +15% / -10% (207 V ... 264 V)

line frequency 50 ... 60 Hz.

The line voltage setting can be changed from the rear without opening the unit.



1. Turn off the controller and disconnect the mains cable.

2. The fuse holder (R10) is located below the 3-pole power connector of the mains jack (R9). Release the fuse holder by pressing its plastic retainers with the aid of a small screwdriver. The retainers are located on the right and left side of the holder and must be pressed towards the center.

3. Unplug the white line voltage switch/indicator (R8, containing the left fuse) from the fuse holder (R10), rotate it until the appropriate voltage marking (100V, 115V, or 230V) is on target for the cutout (R12) of the fuse holder, and plug it back into the fuse holder. Press in the fuse holder until locked on both sides. The appropriate line

voltage marking must be visible in the cutout (R12) of the fuse holder.

### Attention


If you have changed to or from 230 V, also change the mains fuses to the correct value given in section [Replacing Mains Fuses](#).

## 4.2 Replacing Mains Fuses

The two power input fuses are externally accessible. If they blew due to line distortions, incorrect line voltage or other causes, they can be replaced from the rear without opening the unit.

### Attention

To avoid risk of fire only the appropriate fuses for the corresponding line voltage must be used



1. Turn off the LDC200C Series and disconnect the mains cable.

2. The fuse holder (R10) is located below the 3-pole power connector of the mains jack (R9). Release the fuse holder by pressing its plastic retainers using a small screwdriver. The retainers are located on the right and left side of the holder and must be pressed towards the center.

3. Replace the defective fuses (R11) and press in the fuse holder until locked on both sides. Take care to maintain the correct rotation of the white line voltage indicator / switch (R8) which contains the left fuse and is plugged into the fuse holder. The appropriate line voltage marking must be visible in

the cutout (R12) of the fuse holder.

### Fuse types

LDC200VC, LDC201CU, LDC202C, LDC205C, LDC210C, and LDC220C:

100 V	500 mA, time-lag, 250V	T0.5A250V
115 V	500 mA, time-lag, 250V	T0.5A250V
230 V	250 mA, time-lag, 250V	T0.25A250V

LDC240C:

100 V	800 mA, time-lag, 250V	T0.8A250V
115 V	800 mA, time-lag, 250V	T0.8A250V
230 V	400 mA, time-lag, 250V	T0.4A250V

All fuses must meet IEC specification 60127-2/III, time characteristic: time-lag (T), 250V AC, size 5 x 20 mm.

## 4.3 Troubleshooting

In case that your LDC200C Series Controller shows malfunction please check the following items:

- ❑ Unit does not work at all (no display at the front):
  - LDC2xxC connected properly to the mains?
    - Check the power cord and the [line voltage setting](#)
  - LDC2xxC turned on?
    - Turn on your LDC200C Series Controller (button "LINE" F10).
  - Check the [fuses](#) at the rear panel.
    - If blown, replace the fuses with the correct type.
- ❑ The desired laser output power cannot be achieved
  - Is the interlock closed?
    - Check the connection of [interlock](#) pins of the connector jack.
  - Is the laser output turned on ("LASER ON" button)?
    - Change the setting from "OFF" to "ON". The LED "ON" on the front panel of the main-frame must light up.
  - Is the current limit  $I_{LM}$  set to 0?
    - Adjust the [hardware limit  \$I\_{LM}\$](#)  (potentiometer on the front panel) to an appropriate value.
  - Is the laser diode installed properly?
    - Check the connection cable.
  - Is the polarity of the laser diode set correctly?
    - If not change the polarity of the diode or select on the LDC2xxC the [opposite polarity](#).
  - Is the photodiode connected properly?
    - Check the connection cable.
  - Is the photodiode poled correctly?
    - If in constant power mode the displayed photodiode current is negative, the polarity of the photodiode must be reverted.
  - Are you using a [bias voltage](#) with the photodiode in photo current mode?
    - Change the polarity of the diode for photo element mode or change the polarity of the bias voltage source.
  - Is the desired output power set correctly?
    - Check the optical [power display calibration](#)
    - Adjust the desired output power  $P_{LD}$  using the tuning knob.



- The unit switches on, but display shows error message (e.g., “Err06”)
  - This indicates a malfunction of the LDC200C Series Controller.
    - In such case, the controller needs to be returned to Thorlabs for maintenance. Please contact Thorlabs with the information of the error code number and the serial number of your LDC200C Series Controller in order to receive the RMA (Return Material Authorization) instructions accordingly.

If you don't find the error source by means of the trouble shooting list please contact Thorlabs for advise and/or return instructions

## 5 Appendix

### 5.1 Technical Data

#### Common Technical Data

General Data	
Safety Features	Interlock Laser Current Limit Soft Start Short Circuit when Laser off Open Circuit Detection Overtemperature Protection
Display	LED, 5 Digits
Connector for Laser, Photodiode, Interlock & Laser On Signal	9-pin D-Sub Jack
Connectors for Control Input / Output	BNC
Chassis Ground Connector	4mm Banana Jack
Line Voltage	100 V 115 V +15% / -10% 230 V
Mains Supply Overvoltage	Category II (Cat II)
Line Frequency	50 to 60 Hz
Operating Temperature Range <sup>1)</sup>	0°C to +40 °C
Storage Temperature Range	-40°C to +70 °C
Relative Humidity	Max. 80% up to 31 °C, decreasing to 50% at 40 °C
Pollution Degree (Indoor Use only)	2
Operation Altitude	<2000 m
Warm-up Time for Rated Accuracy	10 min
Dimensions (W x H x D) w/o Operating Elements with Operating Elements	146 x 66 x 290 mm <sup>3</sup> 146 x 77 x 320 mm <sup>3</sup>

<sup>1)</sup> non-condensing

All technical data are valid at 23 ± 5°C and 45 ± 15% rel. humidity (non condensing)

## 5.1.1 LDC200CV

<b>Current Control (Constant Current Mode)</b>	
Control Range	0 to $\pm 20$ mA
Compliance Voltage	>6 V
Resolution	1.0 $\mu$ A
Accuracy	$\pm 20$ $\mu$ A
Noise Without Ripple (10Hz to 10MHz, rms, typ.)	<1.0 $\mu$ A
Ripple (50/60Hz, rms, typ.)	<0.5 $\mu$ A
Transients (typ.)	<10 $\mu$ A
Drift, 24hours (typ., 0-10Hz, at constant ambient temperature)	<1 $\mu$ A
Temperature Coefficient	<50 ppm/ $^{\circ}$ C
<b>Current Limit</b>	
Setting Range	0 to >20 mA
Resolution	1 $\mu$ A
Accuracy	$\pm 50$ $\mu$ A
<b>Power Control (Constant Power Mode)</b>	
Photo Current Control Range	5 $\mu$ A to 2 mA
Photo Current Resolution	0.1 $\mu$ A
Photo Current Accuracy	$\pm 2$ $\mu$ A
<b>Analog Modulation Input</b>	
Input Resistance	10 k $\Omega$
Small Signal 3dB Bandwidth, CC Mode	DC to 100 kHz
Modulation Coefficient, CC Mode	2 mA/V $\pm 5\%$
Modulation Coefficient, CP Mode	0.2 mA/V $\pm 5\%$
<b>Laser Current Monitor Output</b>	
Load Resistance	>10 k $\Omega$
Transmission Coefficient	500 V/A $\pm 5\%$
<b>General data</b>	
Maximum Power Consumption	20 VA
Weight	<3.1 kg

All technical data are valid at  $23 \pm 5^{\circ}$ C and  $45 \pm 15\%$  humidity

## 5.1.2 LDC201CU

<b>Current Control (Constant Current Mode)</b>	
Control Range	0 to $\pm 100$ mA
Compliance Voltage	>5 V
Resolution	10 $\mu$ A
Accuracy	$\pm 50$ $\mu$ A
Noise Without Ripple (10Hz to 10MHz, rms, typ.)	<0.2 $\mu$ A
Ripple (50/60Hz, rms, typ.)	<0.5 $\mu$ A
Transients (typ.)	<10 $\mu$ A
Drift, 24hours (typ., 0-10Hz, at constant ambient temperature)	<2 $\mu$ A
Temperature Coefficient	<50 ppm/ $^{\circ}$ C
<b>Current Limit</b>	
Setting Range	0 to >100 mA
Resolution	10 $\mu$ A
Accuracy	$\pm 200$ $\mu$ A
<b>Power Control (Constant Power Mode)</b>	
Photo Current Control Range	25 $\mu$ A to 10 mA
Photo Current Resolution	1 $\mu$ A
Photo Current Accuracy	$\pm 10$ $\mu$ A
<b>Analog Modulation Input</b>	
Input Resistance	10 k $\Omega$
Small Signal 3dB Bandwidth, CC Mode	DC to 0.2 kHz
Modulation Coefficient, CC Mode	10 mA/V $\pm 5\%$
Modulation Coefficient, CP Mode	1 mA/V $\pm 5\%$
<b>Laser Current Monitor Output</b>	
Load Resistance	>10 k $\Omega$
Transmission Coefficient	100 V/A $\pm 5\%$
<b>General data</b>	
Maximum Power Consumption	20 VA
Weight	<3.1 kg

All technical data are valid at  $23 \pm 5^{\circ}$ C and  $45 \pm 15\%$  humidity

### 5.1.3 LDC202C

<b>Current Control (Constant Current Mode)</b>	
Control Range	0 to $\pm 200$ mA
Compliance Voltage	$>10$ V
Resolution	$10$ $\mu$ A
Accuracy	$\pm 100$ $\mu$ A
Noise Without Ripple (10Hz to 10MHz, rms, typ.)	$<1.5$ $\mu$ A
Ripple (50/60Hz, rms, typ.)	$<1.5$ $\mu$ A
Transients (typ.)	$<0.2$ mA
Drift, 24hours (typ., 0-10Hz, at constant ambient temperature)	$<3$ $\mu$ A
Temperature Coefficient	$<50$ ppm/ $^{\circ}$ C
<b>Current Limit</b>	
Setting Range	0 to $>200$ mA
Resolution	$10$ $\mu$ A
Accuracy	$\pm 500$ $\mu$ A
<b>Power Control (Constant Power Mode)</b>	
Photo Current Control Range	$25$ $\mu$ A to $10$ mA
Photo Current Resolution	$1$ $\mu$ A
Photo Current Accuracy	$\pm 10$ $\mu$ A
<b>Analog Modulation Input</b>	
Input Resistance	$10$ k $\Omega$
Small Signal 3dB Bandwidth, CC Mode	DC to $250$ kHz
Modulation Coefficient, CC Mode	$20$ mA/V $\pm 5\%$
Modulation Coefficient, CP Mode	$1$ mA/V $\pm 5\%$
<b>Laser Current Monitor Output</b>	
Load Resistance	$>10$ k $\Omega$
Transmission Coefficient	$50$ V/A $\pm 5\%$
<b>General data</b>	
Maximum Power Consumption	$25$ VA
Weight	$<3.1$ kg

All technical data are valid at  $23 \pm 5^{\circ}$ C and  $45 \pm 15\%$  humidity

## 5.1.4 LDC205C

<b>Current Control (Constant Current Mode)</b>	
Control Range	0 to $\pm 500$ mA
Compliance Voltage	$>10$ V
Resolution	$10$ $\mu$ A
Accuracy	$\pm 0.5$ mA
Noise Without Ripple (10Hz to 10MHz, rms, typ.)	$< 3$ $\mu$ A
Ripple (50/60Hz, rms, typ.)	$< 2$ $\mu$ A
Transients (typ.)	$< 0.5$ mA
Drift, 24hours (typ., 0-10Hz, at constant ambient temperature)	$<10$ $\mu$ A
Temperature Coefficient	$<50$ ppm/ $^{\circ}$ C
<b>Current Limit</b>	
Setting Range	0 to $>500$ mA
Resolution	$10$ $\mu$ A
Accuracy	$\pm 1.5$ mA
<b>Power Control (Constant Power Mode)</b>	
Photo Current Control Range	$25$ $\mu$ A to $10$ mA
Photo Current Resolution	$1$ $\mu$ A
Photo Current Accuracy	$\pm 10$ $\mu$ A
<b>Analog Modulation Input</b>	
Input Resistance	$10$ k $\Omega$
Small Signal 3dB Bandwidth, CC Mode	DC to $150$ kHz
Modulation Coefficient, CC Mode	$50$ mA/V $\pm 5\%$
Modulation Coefficient, CP Mode	$1$ mA/V $\pm 5\%$
<b>Laser Current Monitor Output</b>	
Load Resistance	$>10$ k $\Omega$
Transmission Coefficient	$20$ V/A $\pm 5\%$
<b>General data</b>	
Maximum Power Consumption	$30$ VA
Weight	$<3.1$ kg

All technical data are valid at  $23 \pm 5^{\circ}$ C and  $45 \pm 15\%$  humidity

## 5.1.5 LDC210C

<b>Current Control (Constant Current Mode)</b>	
Control Range	0 to $\pm 1$ A
Compliance Voltage	>10 V
Resolution	100 $\mu$ A
Accuracy	$\pm 1.0$ mA
Noise Without Ripple (10Hz to 10MHz, rms, typ.)	<5 $\mu$ A
Ripple (50/60Hz, rms, typ.)	<3 $\mu$ A
Transients (typ.)	<1 mA
Drift, 24hours (typ., 0-10Hz, at constant ambient temperature)	<20 $\mu$ A
Temperature Coefficient	<50 ppm/ $^{\circ}$ C
<b>Current Limit</b>	
Setting Range	0 to >1 A
Resolution	100 $\mu$ A
Accuracy	$\pm 2.5$ mA
<b>Power Control (Constant Power Mode)</b>	
Photo Current Control Range	50 $\mu$ A to 20 mA
Photo Current Resolution	1 $\mu$ A
Photo Current Accuracy	$\pm 20$ $\mu$ A
<b>Analog Modulation Input</b>	
Input Resistance	10 k $\Omega$
Small Signal 3dB Bandwidth, CC Mode	DC to 100 kHz
Modulation Coefficient, CC Mode	100 mA/V $\pm 5\%$
Modulation Coefficient, CP Mode	2 mA/V $\pm 5\%$
<b>Laser Current Monitor Output</b>	
Load Resistance	>10 k $\Omega$
Transmission Coefficient	10 V/A $\pm 5\%$
<b>General data</b>	
Maximum Power Consumption	40 VA
Weight	<3.1 kg

All technical data are valid at  $23 \pm 5^{\circ}$ C and  $45 \pm 15\%$  humidity

## 5.1.6 LDC220C

<b>Current Control (Constant Current Mode)</b>	
Control Range	0 to $\pm 2$ A
Compliance Voltage	>4 V
Resolution	100 $\mu$ A
Accuracy	$\pm 2.0$ mA
Noise Without Ripple (10Hz to 10MHz, rms, typ.)	<15 $\mu$ A
Ripple (50/60Hz, rms, typ.)	<5 $\mu$ A
Transients (typ.)	<2 mA
Drift, 24hours (typ., 0-10Hz, at constant ambient temperature)	<100 $\mu$ A
Temperature Coefficient	<50 ppm/ $^{\circ}$ C
<b>Current Limit</b>	
Setting Range	0 to >2 A
Resolution	100 $\mu$ A
Accuracy	$\pm 5$ mA
<b>Power Control (Constant Power Mode)</b>	
Photo Current Control Range	50 $\mu$ A to 20 mA
Photo Current Resolution	1 $\mu$ A
Photo Current Accuracy	$\pm 20$ $\mu$ A
<b>Analog Modulation Input</b>	
Input Resistance	10 k $\Omega$
Small Signal 3dB Bandwidth, CC Mode	DC to 50 kHz
Modulation Coefficient, CC Mode	200 mA/V $\pm 5\%$
Modulation Coefficient, CP Mode	2 mA/V $\pm 5\%$
<b>Laser Current Monitor Output</b>	
Load Resistance	>10 k $\Omega$
Transmission Coefficient	5 V/A $\pm 5\%$
<b>General data</b>	
Maximum Power Consumption	60 VA
Weight	<3.3 kg

All technical data are valid at  $23 \pm 5^{\circ}$ C and  $45 \pm 15\%$  humidity



## 5.1.7 LDC240C

<b>Current Control (Constant Current Mode)</b>	
Control Range	0 to $\pm 4$ A
Compliance Voltage	$>5$ V
Resolution	100 $\mu$ A
Accuracy	$\pm 4.0$ mA
Noise Without Ripple (10Hz to 10MHz, rms, typ.)	$<50$ $\mu$ A
Ripple (50/60Hz, rms, typ.)	$<8$ $\mu$ A
Transients (typ.)	$<4$ mA
Drift, 24hours (typ., 0-10Hz, at constant ambient temperature)	$<200$ $\mu$ A
Temperature Coefficient	$<50$ ppm/ $^{\circ}$ C
<b>Current Limit</b>	
Setting Range	0 to $>4$ A
Resolution	100 $\mu$ A
Accuracy	$\pm 10$ mA
<b>Power Control (Constant Power Mode)</b>	
Photo Current Control Range	50 $\mu$ A to 20 mA
Photo Current Resolution	1 $\mu$ A
Photo Current Accuracy	$\pm 20$ $\mu$ A
<b>Analog Modulation Input</b>	
Input Resistance	10 k $\Omega$
Small Signal 3dB Bandwidth, CC Mode	DC to 30 kHz
Modulation Coefficient, CC Mode	400 mA/V $\pm 5\%$
Modulation Coefficient, CP Mode	2 mA/V $\pm 5\%$
<b>Laser Current Monitor Output</b>	
Load Resistance	$>10$ k $\Omega$
Transmission Coefficient	2.5 V/A $\pm 5\%$
<b>General data</b>	
Maximum Power Consumption	100 VA
Weight	$<3.3$ kg

All technical data are valid at  $23 \pm 5^{\circ}$ C and  $45 \pm 15\%$  humidity

## 5.2 Declaration of Conformity

### *EU Declaration of Conformity*

*in accordance with EN ISO 17050-1:2010*

We: Thorlabs GmbH

Of: Hans-Boeckler-Str. 6, 85221 Dachau/München, Deutschland

*in accordance with the following Directive(s):*

2014/35/EU	Low Voltage Directive (LVD)
2014/30/EU	Electromagnetic Compatibility (EMC) Directive
2011/65/EU	Restriction of Use of Certain Hazardous Substances (RoHS)

*hereby declare that:*

Model: **LDC2xxCx**

Equipment: **Benchtop Laserdiode Current Controller**

*is in conformity with the applicable requirements of the following documents:*

EN 61010-1	Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use.	2010
EN 61326-1	Electrical Equipment for Measurement, Control and Laboratory Use - EMC Requirements	2013

*and which, issued under the sole responsibility of Thorlabs, is in conformity with Directive 2011/65/EU of the European Parliament and of the Council of 8th June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, for the reason stated below:*

does not contain substances in excess of the maximum concentration values tolerated by weight in homogenous materials as listed in Annex II of the Directive

*I hereby declare that the equipment named has been designed to comply with the relevant sections of the above referenced specifications, and complies with all applicable Essential Requirements of the Directives.*

Signed:  On: 26 June 2015

Name: Dorothee Jennrich

Position: General Manager

EDC - LDC2xxCx -2015-06-26

**CE 15**

This product was tested for and complies with the following standards:

- CAN/CSA-C22.2 No. 61010-1-04
- ANSI/UL 61010-1-2004

## 5.3 Warranty

Thorlabs warrants material and production of the LDC200C Series for a period of 24 months starting with the date of shipment. During this warranty period Thorlabs will see to defaults by repair or by exchange if these are entitled to warranty.

For warranty repairs or service the unit must be sent back to Thorlabs. The customer will carry the shipping costs to Thorlabs, in case of warranty repairs Thorlabs will carry the shipping costs back to the customer.

If no warranty repair is applicable the customer also has to carry the costs for back shipment.

In case of shipment from outside EU duties, taxes etc. which should arise have to be carried by the customer.

Thorlabs warrants the hard- and/or software determined by Thorlabs for this unit to operate fault-free provided that they are handled according to our requirements. However, Thorlabs does not warrant a fault free and uninterrupted operation of the unit, of the software or firmware for special applications nor this instruction manual to be error free. Thorlabs is not liable for consequential damages.

### Restriction of Warranty

The warranty mentioned before does not cover errors and defects being the result of improper treatment, software or interface not supplied by us, modification, misuse or operation outside the defined ambient stated by us or unauthorized maintenance.

Further claims will not be consented to and will not be acknowledged. Thorlabs does explicitly not warrant the usability or the economical use for certain cases of application.

Thorlabs reserves the right to change this instruction manual or the technical data of the described unit at any time.

## 5.4 Copyright and Exclusion of Reliability

*Thorlabs* has taken every possible care in preparing this document. We however assume no liability for the content, completeness or quality of the information contained therein. The content of this document is regularly updated and adapted to reflect the current status of the hardware and/or software. We furthermore do not guarantee that this product will function without errors, even if the stated specifications are adhered to.

Under no circumstances can we guarantee that a particular objective can be achieved with the purchase of this product.

Insofar as permitted under statutory regulations, we assume no liability for direct damage, indirect damage or damages suffered by third parties resulting from the purchase of this product. In no event shall any liability exceed the purchase price of the product.

Please note that the content of this document is neither part of any previous or existing agreement, promise, representation or legal relationship, nor an alteration or amendment thereof. All obligations of *Thorlabs* result from the respective contract of sale, which also includes the complete and exclusively applicable warranty regulations. These contractual warranty regulations are neither extended nor limited by the information contained in this document. Should you require further information on this product, or encounter specific problems that are not discussed in sufficient detail in the document, please contact your local *Thorlabs* dealer or system installer.

All rights reserved. This document may not be reproduced, transmitted or translated to another language, either as a whole or in parts, without the prior written permission of *Thorlabs*.

Copyright © Thorlabs 2018. All rights reserved.

## 5.5 Thorlabs 'End of Life' Policy (WEEE)

As required by the WEEE (Waste Electrical and Electronic Equipment Directive) of the European Community and the corresponding national laws, Thorlabs offers all end users in the EC the possibility to return “end of life” units without incurring disposal charges.

This offer is valid for Thorlabs electrical and electronic equipment

- sold after August 13<sup>th</sup> 2005
- marked correspondingly with the crossed out “wheelie bin” logo (see figure below)
- sold to a company or institute within the EC
- currently owned by a company or institute within the EC
- still complete, not disassembled and not contaminated

As the WEEE directive applies to self contained operational electrical and electronic products, this “end of life” take back service does not refer to other Thorlabs products, such as

- pure OEM products, that means assemblies to be built into a unit by the user (e. g. OEM laser driver cards)
- components
- mechanics and optics
- left over parts of units disassembled by the user (PCB's, housings etc.).

### Waste treatment on your own responsibility

If you do not return an “end of life” unit to Thorlabs, you must hand it to a company specialized in waste recovery. Do not dispose of the unit in a litter bin or at a public waste disposal site.

WEEE Number (Germany) : DE97581288

### Ecological background

It is well known that waste treatment pollutes the environment by releasing toxic products during decomposition. The aim of the European RoHS Directive is to reduce the content of toxic substances in electronic products in the future.

The intent of the WEEE Directive is to enforce the recycling of WEEE. A controlled recycling of end-of-life products will thereby avoid negative impacts on the environment.



*Crossed out  
"Wheelie Bin" symbol*

## 5.6 Thorlabs Worldwide Contacts

For technical support or sales inquiries, please visit us at [www.thorlabs.com/contact](http://www.thorlabs.com/contact) for our most up-to-date contact information.



### **USA, Canada, and South America**

Thorlabs, Inc.  
[sales@thorlabs.com](mailto:sales@thorlabs.com)  
[techsupport@thorlabs.com](mailto:techsupport@thorlabs.com)

### **Europe**

Thorlabs GmbH  
[europe@thorlabs.com](mailto:europe@thorlabs.com)

### **France**

Thorlabs SAS  
[sales.fr@thorlabs.com](mailto:sales.fr@thorlabs.com)

### **Japan**

Thorlabs Japan, Inc.  
[sales@thorlabs.jp](mailto:sales@thorlabs.jp)

### **UK and Ireland**

Thorlabs Ltd.  
[sales.uk@thorlabs.com](mailto:sales.uk@thorlabs.com)  
[techsupport.uk@thorlabs.com](mailto:techsupport.uk@thorlabs.com)

### **Scandinavia**

Thorlabs Sweden AB  
[scandinavia@thorlabs.com](mailto:scandinavia@thorlabs.com)

### **Brazil**

Thorlabs Vendas de Fotônicos Ltda.  
[brasil@thorlabs.com](mailto:brasil@thorlabs.com)

### **China**

Thorlabs China  
[chinasales@thorlabs.com](mailto:chinasales@thorlabs.com)





**THORLABS**

[www.thorlabs.com](http://www.thorlabs.com)

---